

DOC NUMBER	TITLE	REFERENCE SOP
PHDB1/BIO/FOR-28.5	DAILY AND WEEKLY MAINTENANCE FOR COBAS 6000	PHDB1/BIO/SOP-028
ISSUED DATE: 01-07-2018		Version 1.0
EFFECTIVE DATE: 01-07-2018		Supersedes version no.: N/A

Equipment number: **PHD1-BIO-MAA01-029**

Month: _____ Year: _____

DAILY CHECKLIST	01	02	03	04	05	06	07	08	09	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Clean Instrument Surfaces *																															
Check Waste Solution reservoir																															
Check Paper Supply																															
C501 MODULE																															
Process Green Rack																															
Clean sample probe & shield																															
Pipe**																															
ISE sipper probe**																															
Reagent Probes**																															
Cell Rinse nozzles***																															
Drain port ***																															
Perform Daily Pipe																															
Photometer check reading at 340 nm (Reading < 14000)																															
WEEKLY CHECKLIST																															
System Power Off																															
System Power On																															
Perform weekly pipe																															
Clean Cell Covers**																															
Clean Rinse Stations****																															
Clean IS bath***																															
TECH INITIALS																															

Location: Biochemistry

Page 1 of 2

Proprietary information: This document and the contents hereof are considered proprietary and confidential information of PHDB1 Laboratory and disclosure to unauthorized individuals or dissemination, publication or copying is prohibited without prior written consent.

CONTROLLED DOCUMENT:

DOC NUMBER	TITLE	REFERENCE SOP
PHDB1/BIO/FOR-28.5	DAILY AND WEEKLY MAINTENANCE FOR COBAS 6000	PHDB1/BIO/SOP-028
ISSUED DATE: 01-07-2018		Version 1.0
EFFECTIVE DATE: 01-07-2018		Supersedes version no.: N/A

Equipment number: PHD1-BIO-MAA01-029

Month: _____ Year: _____

DAILY CHECKLIST	01	02	03	04	05	06	07	08	09	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
C601 MODULE																															
Clean : Sample probe ***																															
Reagent Probe**																															
Sipper Probes**																															
Pre wash sipper probe**																															
WEEKLY/AS NEEDED																															
Perform Liquid Flow Cleaning																															
Clean PC/CC nozzles, tubes filters and bottle stand***																															
Clean Pre Clean mixer and separation station***																															
Clean Incubator***																															
Clean Assay Cup vortex mixer***																															
Clean microbead mixer*****																															
Clean rinse Station*****																															
Clean reagent disk***																															
Clean solid waste compartments																															
TECH INITIALS																															

*Clean with disinfectant (Cavicide) **Clean with alcohol. *** Clean with Distilled or deionized water ****Clean with 2% ECOtergent followed by water

*****Clean with alcohol followed by water *****Clean with alcohol then 50mL 2% Ecotergent followed by 50 mL DI water

Remarks:

Reviewed by:

Laboratory Director/Designee: _____ Date: _____

Location: Biochemistry

Page 2 of 2